

# BIOMEX

We serve different

Your specialist for biological and clinical raw materials for more than 25 years



Several IvD manufactureres are selling SARS-CoV-2 antibody assays without offering the corresponding internal quality controls. This is mainly due to the lack of available positive plasma on the market. The controls are important to assess laboratory precision and day-to-day assay performance of the testing systems.

To fill this gap, *DiaMex GmbH* (a sister company of *Biomex GmbH*) has developed an Anti-SARS-CoV-2 IgG and IgM control to be used with several different assays. The controls

- are human plasma-based to mimic patient specimen encountered in a daily testing environment and provide a complete quality solution for SARS-CoV-2 antibody testing
- react at an appropriate level within the analytical range of the assays with a high lot to lot consistency
- can be used as third party controls to verify systematic errors and satisfy regulatory guidelines for use of independent controls.



## Ordering information:

DiaMex GmbH  
Siemensstraße 38  
69123 Heidelberg  
[order@diamex.com](mailto:order@diamex.com)

Product code	QC Product	Kit Format	Vial Volume [mL]
SR14039	Optitrol Covid-19 IgG	1 x pos	1 mL
SR14037	Optitrol Covid-19 IgG	2 x pos.	2,5 mL
SR14033	Optitrol Covid-19 IgG	4 x pos.	2,5 mL
SR14049	Optitrol Covid-19 IgM	1 x pos	1 mL
SR14047	Optitrol Covid-19 IgM	2 x pos.	2,5 mL
SR14043	Optitrol Covid-19 IgM	4 x pos.	2,5 mL

# DiaMex

GmbH

Highend diagnostic products

For further information please contact the staff at DiaMex ([order@diamex.com](mailto:order@diamex.com)). Our entire team is here to assist you.